

# Human C1orf56 / Uncharacterized protein C1orf56 ELISA Kit



Catalog No: #EK21050

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human C1orf56 / Uncharacterized protein C1orf56 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	UNQ547/PRO1104
Accession No.	Uniprot:Q9BUN1
Uniprot	Q9BUN1
Storage	6 months

## Application Details

Detection Range:15.6-1000pg/mL Assay Type:Sandwich Sensitivity:8.3pg/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only