

Human TM4SF5 / Transmembrane 4 L6 family member 5 ELISA Kit



Catalog No: #EK21051

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human TM4SF5 / Transmembrane 4 L6 family member 5 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Tetraspan transmembrane protein L6H
Accession No.	Uniprot:O14894
Uniprot	O14894
Storage	6 months

Application Details

Detection Range:0.156-10ng/mL Assay Type:Competitive Sensitivity:0.057ng/mL

Product Description

Detection Method:Competitive

Note: This product is for in vitro research use only