

Mouse Cklf / Chemokine-like factor ELISA Kit

Catalog No: #EK21052



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cklf / Chemokine-like factor ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Accession No.	Uniprot:Q9DAS1
Uniprot	Q9DAS1
Storage	6 months

Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:39.1pg/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only