

Human C10orf91 / Uncharacterized protein C10orf91 ELISA Kit



Catalog No: #EK21053

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human C10orf91 / Uncharacterized protein C10orf91 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Accession No.	Uniprot:Q5T1B1
Uniprot	Q5T1B1
Storage	6 months

Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.12ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only