

Human NR2E3 / Photoreceptor-specific nuclear receptor ELISA Kit

Catalog No: #EK21061



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human NR2E3 / Photoreceptor-specific nuclear receptor ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Nuclear receptor subfamily 2 group E member 3, Retina-specific nuclear receptor, PNR, RNR
Accession No.	Uniprot:Q9Y5X4
Uniprot	Q9Y5X4
Storage	6 months

Application Details

Detection Range: 0.78-50ng/mL Assay Type: Sandwich Sensitivity: 0.395ng/mL

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only