

# Human FAF1 / FAS-associated factor 1 ELISA Kit

Catalog No: #EK21065



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human FAF1 / FAS-associated factor 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	hFAF1,UBX domain-containing protein 12,UBX domain-containing protein 3A,UBXD12,UBXN3A,CGI-03
Accession No.	Uniprot:Q9UNN5
Uniprot	Q9UNN5
Storage	6 months

## Application Details

Detection Range:0.78-50ng/mL Assay Type:Sandwich Sensitivity:0.375ng/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only