

Human FOLH1 / Glutamate carboxypeptidase 2 ELISA Kit

Catalog No: #EK21067

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human FOLH1 / Glutamate carboxypeptidase 2 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Cell growth-inhibiting gene 27 protein,Folate hydrolase 1,Folypoly-gamma-glutamate carboxypeptidase,FGCP,Glutamate carboxypeptidase II,GCPII,Membrane glutamate carboxypeptidase,mGCP,N-acetylated-alpha-linked acidic dipeptidase I,NAALADase I,Prostate-specific membrane antigen,PSM,Pteroylpoly-gamma-glutamate carboxypeptidase,FOLH,NAALAD1,PSM,PSMA,GIG27
Accession No.	Uniprot:Q04609
Uniprot	Q04609
Storage	6 months

Application Details

Detection Range:0.625-40ng/mL Assay Type:Sandwich Sensitivity:0.295ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only