

# Human ZC3H8 / Zinc finger CCCH domain-containing protein 8 ELISA Kit



Catalog No: #EK21068

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ZC3H8 / Zinc finger CCCH domain-containing protein 8 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	ZC3HDC8
Accession No.	Uniprot:Q8N5P1
Uniprot	Q8N5P1
Storage	6 months

## Application Details

Detection Range:0.156-10ng/mL Assay Type:Sandwich Sensitivity:0.057ng/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only