

# Human CCAR2 / Cell cycle and apoptosis regulator protein 2 ELISA Kit



Catalog No: #EK21072

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Human CCAR2 / Cell cycle and apoptosis regulator protein 2 ELISA Kit  |
| Brief Description  | ELISA Kit   |
| Applications       | ELISA   |
| Species Reactivity | Human   |
| Other Names        | Cell division cycle and apoptosis regulator protein 2, DBIRD complex subunit KIAA1967, Deleted in breast cancer gene 1 protein, DBC-1, NET35, p30 DBC, DBC1, KIAA1967 |
| Accession No.      | Uniprot:Q8N163  |
| Uniprot            | Q8N163  |
| Storage            | 6 months  |

## Application Details

Detection Range: 0.78-50ng/mL Assay Type: Sandwich Sensitivity: 0.33ng/mL

## Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only