

# Human TBXA2R / Thromboxane A2 receptor ELISA Kit



Catalog No: #EK21073

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human TBXA2R / Thromboxane A2 receptor ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	TXA2-R, Prostanoid TP receptor
Accession No.	Uniprot:P21731
Uniprot	P21731
Storage	6 months

## Application Details

Detection Range: 0.156-10ng/mL Assay Type: Sandwich Sensitivity: 0.039ng/mL

## Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only