

# Human MT-ND1 / NADH-ubiquinone oxidoreductase chain 1 ELISA Kit



Catalog No: #EK21074

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Human MT-ND1 / NADH-ubiquinone oxidoreductase chain 1 ELISA Kit |
| Brief Description  | ELISA Kit   |
| Applications       | ELISA   |
| Species Reactivity | Human   |
| Other Names        | NADH dehydrogenase subunit 1,MTND1,NADH1,ND1                    |
| Accession No.      | Uniprot:P03886  |
| Uniprot            | P03886  |
| Storage            | 6 months  |

## Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.144ng/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only