

Human FKRP / Fukutin-related protein ELISA Kit

Catalog No: #EK21075



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human FKRP / Fukutin-related protein ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Ribitol-5-phosphate transferase
Accession No.	Uniprot:Q9H9S5
Uniprot	Q9H9S5
Storage	6 months

Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.13ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only