

# Human ALAS2 / 5-aminolevulinate synthase, erythroid-specific, mitochondrial ELISA Kit



Catalog No: #EK21077

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ALAS2 / 5-aminolevulinate synthase, erythroid-specific, mitochondrial ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	ALAS-E,5-aminolevulinic acid synthase 2,Delta-ALA synthase 2,Delta-aminolevulinate synthase 2,ALASE,ASB
Accession No.	Uniprot:P22557
Uniprot	P22557
Storage	6 months

## Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.12ng/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only