

Human SH3RF1 / E3 ubiquitin-protein ligase SH3RF1 ELISA Kit



Catalog No: #EK21079

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human SH3RF1 / E3 ubiquitin-protein ligase SH3RF1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Plenty of SH3s,Protein POSH,RING finger protein 142,RING-type E3 ubiquitin transferase SH3RF1,SH3 domain-containing RING finger protein 1,SH3 multiple domains protein 2,KIAA1494,POSH,POSH1,RNF142,SH3MD2
Accession No.	Uniprot:Q7Z6J0
Uniprot	Q7Z6J0
Storage	6 months

Application Details

Detection Range:0.78-50ng/mL Assay Type:Sandwich Sensitivity:0.365ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only