

Human ECD / Protein ecdysoneless homolog ELISA Kit



Catalog No: #EK21080

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human ECD / Protein ecdysoneless homolog ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Human suppressor of GCR two,hSGT1
Accession No.	Uniprot:O95905
Uniprot	O95905
Storage	6 months

Application Details

Detection Range:31.2-2000pg/mL Assay Type:Sandwich Sensitivity:12.5 pg/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only