

Human SNRPB / Small nuclear ribonucleoprotein-associated proteins B and B' ELISA Kit

Catalog No: #EK21082

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human SNRPB / Small nuclear ribonucleoprotein-associated proteins B and B' ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	snRNP-B, Sm protein B/B', Sm-B/B', COD, SNRPB1
Accession No.	Uniprot:P14678
Uniprot	P14678
Storage	6 months

Application Details

Detection Range:0.156-10ng/mL Assay Type:Sandwich Sensitivity:0.051ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only