

Human BRD2 / Bromodomain-containing protein 2 ELISA Kit



Catalog No: #EK21083

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human BRD2 / Bromodomain-containing protein 2 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Really interesting new gene 3 protein, KIAA9001, RING3, O27.1.1
Accession No.	Uniprot:P25440
Uniprot	P25440
Storage	6 months

Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.11ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only