

Rat Opa1 / Dynamin-like 120 kDa protein, mitochondrial ELISA Kit



Catalog No: #EK21084

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Opa1 / Dynamin-like 120 kDa protein, mitochondrial ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rt
Other Names	Optic atrophy protein 1 homolog
Accession No.	Uniprot:Q2TA68
Uniprot	Q2TA68
Storage	6 months

Application Details

Detection Range:0.156-10ng/mL Assay Type:Sandwich Sensitivity:0.054ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only