

Human PEX1 / Peroxisome biogenesis factor 1 ELISA Kit

Catalog No: #EK21086



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Human PEX1 / Peroxisome biogenesis factor 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Peroxin-1,Peroxisome biogenesis disorder protein 1
Accession No.	Uniprot:O43933
Uniprot	O43933
Storage	6 months

Application Details

Detection Range:0.312-20ng/mLAssay Type:SandwichSensitivity:0.135ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only