

# Human PAWR / PRKC apoptosis WT1 regulator protein ELISA Kit

Catalog No: #EK21087



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human PAWR / PRKC apoptosis WT1 regulator protein ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Prostate apoptosis response 4 protein,Par-4,PAR4
Accession No.	Uniprot:Q96IZ0
Uniprot	Q96IZ0
Storage	6 months

## Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:35pg/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only