

Human TMEM201 / Transmembrane protein 201 ELISA Kit



Catalog No: #EK21089

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human TMEM201 / Transmembrane protein 201 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Spindle-associated membrane protein 1,NET5,SAMP1
Accession No.	Uniprot:Q5SNT2
Uniprot	Q5SNT2
Storage	6 months

Application Details

Detection Range:0.156-10ng/mL Assay Type:Sandwich Sensitivity:0.055ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only