

Human RPGR / X-linked retinitis pigmentosa GTPase regulator ELISA Kit



Catalog No: #EK21090

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human RPGR / X-linked retinitis pigmentosa GTPase regulator ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	RP3,XLRP3
Accession No.	Uniprot:Q92834
Uniprot	Q92834
Storage	6 months

Application Details

Detection Range:0.78-50ng/mL Assay Type:Sandwich Sensitivity:0.33ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only