

## Mouse Nfkbiz / NF-kappa-B inhibitor zeta ELISA Kit

Catalog No: #EK21092

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Nfkbiz / NF-kappa-B inhibitor zeta ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Other Names	I-kappa-B-zeta,IkB-zeta,IL-1 inducible nuclear ankyrin-repeat protein,INAP,Molecule possessing ankyrin repeats induced by lipopolysaccharide,MAIL,Inap,Mail
Accession No.	Uniprot:Q9EST8
Uniprot	Q9EST8
Storage	6 months

## Application Details

Detection Range:0.156-10ng/ml Assay Type:Sandwich Sensitivity:0.059ng/ml

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only