

# Human NPR1 / Atrial natriuretic peptide receptor 1 ELISA Kit



Catalog No: #EK21100

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human NPR1 / Atrial natriuretic peptide receptor 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Atrial natriuretic peptide receptor type A, ANP-A, Guanylate cyclase A, GC-A, ANPRA
Accession No.	Uniprot:P16066
Uniprot	P16066
Storage	6 months

## Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.15ng/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only