

# Human ACE2 / Angiotensin-converting enzyme 2 ELISA Kit



Catalog No: #EK21102

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ACE2 / Angiotensin-converting enzyme 2 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	ACE-related carboxypeptidase,Angiotensin-converting enzyme homolog,ACEH,Metalloprotease MPROT15,UNQ868/PRO1885
Accession No.	Uniprot:Q9BYF1
Uniprot	Q9BYF1
Storage	6 months

## Application Details

Detection Range:0.78-50ng/mL Assay Type:Sandwich Sensitivity:0.34ng/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only