

# Mouse LCMV ELISA Kit

Catalog No: #EK22003



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse LCMV ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

## Application Details

Detect Range: Calculate Critical (CUT OFF): Critical = the average of Negative control well + 0.15. Negative Result: sample OD < Calculate Critical (CUT OFF) is Negative. Positive Result: sample OD > Calculate Critical (CUT OFF) is Positive.

Sample Type: Serum, plasma, Urine, tissue homogenates and other biological fluids.

Sample Volume: 50ul

Assay Time: 3-5h

Detection wavelength: 450 nm

Note: This product is for in vitro research use only