

Human CD155 Protein, mFc-His tag

Catalog No: #AP89554

Package Size: #AP89554-1 10ug #AP89554-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CD155 Protein, mFc-His tag
Host Species	HEK293
Purification	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Species Reactivity	Human
Immunogen Description	CD155(Trp21-Asn343)+mFc(Pro99-Lys330)+6xHis tag
Other Names	PVR, FLJ25946, PVS, CD155, TAGE4, HVED, NECL5
Calculated MW	62.2 kDa
Tag Info	C-Mouse Fc and 6xHis tag
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 34 % trehalose is added as protectants before lyophilization.
Storage	Store at -80°C for 12 months (Avoid repeated freezing and thawing)

Images




Figure 1. Human CD155 Protein, mFc-His Tag on SDS-PAGE under reducing condition.



Product Description

Expression Region:550Research Topic:The protein encoded by this gene is a transmembrane glycoprotein belonging to the immunoglobulin superfamily. The external domain mediates cell attachment to the extracellular matrix molecule vitronectin, while its intracellular domain interacts with

the dynein light chain Tctex-1/DYNLT1. The gene is specific to the primate lineage, and serves as a cellular receptor for poliovirus in the first step of poliovirus replication. Multiple transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only