

Human GPRC5D Protein, hFc-His Tag

Catalog No: #AP89573

Package Size: #AP89573-1 10ug #AP89573-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human GPRC5D Protein, hFc-His Tag
Host Species	HEK293
Purification	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Species Reactivity	Human
Immunogen Description	GPRC5D(Met1-Glu27)+hFc(Glu99-Ala330)+6xHis tag
Other Names	GPRC5D
Calculated MW	31.5 kDa
Tag Info	C-Human Fc and 6xHis tag
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 53% trehalose is added as protectants before lyophilization.
Storage	Store at -80°C for 12 months (Avoid repeated freezing and thawing)

Images





Figure 1. Human GPRC5D Protein, hFc-His Tag on SDS-PAGE under reducing condition.



Product Description

Expression Region:550Research Topic:The protein encoded by this gene is a member of the G protein-coupled receptor family; however, the specific function of this gene has not yet been determined.

Note: This product is for in vitro research use only