

Human TweakR Protein, mFc Tag

Catalog No: #AP89594

Package Size: #AP89594-1 10ug #AP89594-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human TweakR Protein, mFc Tag
Host Species	HEK293
Purification	The purity of the protein is greater than 90% as determined by SDS-PAGE and Coomassie blue staining.
Species Reactivity	Human
Immunogen Description	TweakR(Glu28-Pro80)+mFc(Pro99-Lys330)
Other Names	TNFRSF12A, ?FGF-inducible 14, ?FN14, ?TweakR, ?CD266
Calculated MW	31.8 kDa
Tag Info	C-Mouse Fc tag
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 74 % trehalose is added as protectants before lyophilization.
Storage	Store at -80°C for 12 months (Avoid repeated freezing and thawing)

Images





Figure 1. Human TweakR Protein, mFc Tag on SDS-PAGE under non-reducing(NR) and reducing(R) conditions.



Product Description

Expression Region:510Research Topic:TNFRSF12A (TNF Receptor Superfamily Member 12A) is a Protein Coding gene. Diseases associated with TNFRSF12A include?Multiple Sclerosis?and?Stroke, Ischemic. Among its related pathways are?ERK Signaling?and?CDK-mediated phosphorylation

and removal of Cdc6

Note: This product is for in vitro research use only