

## Human CD34 Protein, His Tag

Catalog No: #AP89599

Package Size: #AP89599-1 10ug #AP89599-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Human CD34 Protein, His Tag
Host Species	HEK293
Purification	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Species Reactivity	Human
Immunogen Description	CD34(Ser32-Thr290)+6xHis tag
Other Names	CD34 Molecule
Calculated MW	28.3 kDa
Tag Info	C- 6xHis tag
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 79 % trehalose is added as protectants before lyophilization.
Storage	Store at -80°C for 12 months (Avoid repeated freezing and thawing)

## Images





Figure 1. Human CD34 Protein, His Tag on SDS-PAGE under reducing condition.



## Product Description

Expression Region:550Research Topic:The protein encoded by this gene may play a role in the attachment of stem cells to the bone marrow extracellular matrix or to stromal cells. This single-pass membrane protein is highly glycosylated and phosphorylated by protein kinase C. Two

transcript variants encoding different isoforms have been found for this gene.

---

Note: This product is for in vitro research use only