

Human GITR Ligand Protein, His tag

Catalog No: #AP89606

Package Size: #AP89606-1 10ug #AP89606-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Human GITR Ligand Protein, His tag |
| Host Species | HEK293 |
| Purification | The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining. |
| Species Reactivity | Human |
| Immunogen Description | GITR Ligand(Gln50-Ser177)+6xHis |
| Other Names | TNFSF18, AITRL, TL6, hGITRL, GITR Ligand |
| Calculated MW | 15.3 kDa |
| Tag Info | C-His tag |
| Formulation | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 86 % trehalose is added as protectants before lyophilization. |
| Storage | Store at -80°C for 12 months (Avoid repeated freezing and thawing) |

Images





Figure 1. Human GITR Ligand Protein, His Tag on SDS-PAGE under reducing condition.



Product Description

Expression Region:600Research Topic:The protein encoded by this gene is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family. This cytokine is a ligand for receptor TNFRSF18/AITR/GITR. It has been shown to modulate T lymphocyte survival in peripheral tissues. This cytokine

is also found to be expressed in endothelial cells, and is thought to be important for interaction between T lymphocytes and endothelial cells.

Note: This product is for in vitro research use only