

Human ICOS Protein, hFc tag

Catalog No: #AP89610

Package Size: #AP89610-1 10ug #AP89610-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human ICOS Protein, hFc tag
Host Species	HEK293
Purification	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Species Reactivity	Human
Immunogen Description	ICOS(Glu21-Lys140)+hFc(Glu99-Ala330)
Other Names	ICOS, CD278, AILIM, Inducible T-cell costimulator
Tag Info	C-Human Fc tag
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 90 % trehalose is added as protectants before lyophilization.
Storage	Store at -80°C for 12 months (Avoid repeated freezing and thawing)

Images





Figure 1. Human ICOS Protein, hFc Tag on SDS-PAGE under reducing condition.



Product Description

Expression Region:600Research Topic:The protein encoded by this gene belongs to the CD28 and CTLA-4 cell-surface receptor family. It forms homodimers and plays an important role in cell-cell signaling, immune responses, and regulation of cell proliferation.

Note: This product is for in vitro research use only