

Human Trop2 Protein, mFc-His Tag

Catalog No: #AP89620

Package Size: #AP89620-1 10ug #AP89620-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Trop2 Protein, mFc-His Tag
Host Species	HEK293
Purification	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Species Reactivity	Human
Immunogen Description	Trop2(His27-Thr274)+mFc(Pro99-Lys330)+6xHis tag
Other Names	TACSTD2,GA733-1,M1S1,TROP2
Calculated MW	54.9 kDa
Tag Info	C-Mouse Fc and 6xHis tag
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 100 % trehalose is added as protectants before lyophilization.
Storage	Store at -80°C for 12 months (Avoid repeated freezing and thawing)

Images





Figure 1. Human Trop2 Protein, mFc-His Tag on SDS-PAGE under reducing condition.



Product Description

Expression Region:550Research Topic:This intronless gene encodes a carcinoma-associated antigen. This antigen is a cell surface receptor that transduces calcium signals. Mutations of this gene have been associated with gelatinous drop-like corneal dystrophy.

Note: This product is for in vitro research use only