

Recombinant human Ras-related protein Rab-3C

Catalog No: #AP71454



Package Size: #AP71454-1 20ug #AP71454-2 100ug #AP71454-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant human Ras-related protein Rab-3C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-227aaSequence Info:Full Length
Accession No.	Q96E17
Uniprot	Q96E17
GeneID	115827;
Calculated MW	42 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	MRHEAPMQMASAQDARYGQKDSSDQNFDFYMFKLLIIGNSSVGKTSFLFRYADDSFTSAFVSTVGIDFKVKTVF KNEKRIKLQIWDTAGQERYRTITTAYYRGAMGFILMYDITNEESFNAVQDWSTQIKTYSWDNAQVILVGNKCD MEDERVISTERGQHLGELGFEFFETSAKDNINVKQTFERLVDIICDKMSESLETDPAITAAKQNTRLKETPPPP QPNCAC
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Protein transport. Probably involved in vesicular traffic .

References

Cloning, mapping, and characterization of the human Rab3C gene.Cheng H., Ma Y., Ni X., Jiang M., Luo Y., Ying K., Xie Y., Ma Y.Biochem. Genet. 40:263-272(2002)Research Topic:Transport

Note: This product is for in vitro research use only