

Recombinant human Paraneoplastic antigen-like protein 5



Catalog No: #AP71481

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP71481-1 20ug #AP71481-2 100ug #AP71481-3 1mg

Description

Product Name	Recombinant human Paraneoplastic antigen-like protein 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:4-442aaSequence Info:Partial
Other Names	Tumor antigen BJ-HCC-25
Accession No.	Q96PV4
Uniprot	Q96PV4
GeneID	114824;
Calculated MW	64.9 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	<p>TLLEDWCKGMDMDPRKALLIVGIPMECSEVEIQDQTVKAGLQPLCAYRVLGRMFRREDNAKAVFIELADTVNYT</p> <p>TLPSHIPGKGGSEWVVKPRNPDDEFLSRLNYFLKDEGRSMTDVARALGCCSLPAESLDAEVMQVRSPPLE</p> <p>PPKESMWYRKLKVFSGTASPSGGEETFEDWLEQVTEIMPIWQVSEVEKRRRLLLESLRGPALSIMRVLQANND</p> <p>SITVEQCLDALKQIFGDKEDFRASQFRFLQTSFKIGEKVSTFLLRLEPLLQAVHKSPLSVRSTDMIRLKHLLAR</p> <p>VAMTPALRGKLELLDQRGCPPNFLELMKLRDEEEWENTEAVMKNKEKPSGRGRGASGRQARAEASVSAPQ</p> <p>ATVQARSFSDSSPQTIQGGPLPVKRRRLLGSESTRGEDHGGATYPKAENQTPGREGPQAAGEELGNEAGA</p> <p>GAMSH</p>
Formulation	Tris-based buffer50% glycerol
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.</p>

References

Cloning and identification of genes which are differentially expressed in carcinoma.Dong X.-Y., Chen W.-F.Prediction of the coding sequences of unidentified human genes. XXI. The complete sequences of 60 new cDNA clones from brain which code for large proteins.Nagase T., Kikuno R., Ohara O.DNA Res. 8:179-187(2001)Research Topic:Cancer

Note: This product is for in vitro research use only