

Recombinant human Endophilin-B2

Catalog No: #AP71542



Package Size: #AP71542-1 20ug #AP71542-2 100ug #AP71542-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant human Endophilin-B2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-395aaSequence Info:Full Length
Other Names	SH3 domain-containing GRB2-like protein B2
Accession No.	Q9NR46
Uniprot	Q9NR46
GeneID	56904;
Calculated MW	60 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	MDFNMKKLASDAGIFFTRAVQFTEEKFGQAEKTELD AHFENLLARADSTKNWTEKILRQTEVLLQPNPSARVE EFLYEKLDKRVPSRVVTNGELLAQYMADAASELGPTTPYGKTLIKVAEAEKQLGAAERDFIHTASISFLTPLRNFL EGDWKTISKERRLLQNRRLDL DACKARLKKAKAAEAKATTVPDFQETRPRNYILSASASALWNDEVDKAEQEL RVAQTEFDRQAEVTRLLLEGISSTHVNHLRCLHEFVKSQTTYAQCYRHMLDLQKQLGRFPGTFVGTTEPASP PLSSTSPPTAAATMPVVPSVASLAPPGEASLCL EEVAPPASGTRKARVLYDYEAADSSSELALLADELITVYSLP GMDPDWLGIRGNKKGKVPVTYLELLS
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

SH3GLB, a new endophilin-related protein family featuring an SH3 domain.Pierrat B., Simonen M., Cueto M., Mestan J., Ferrigno P., Heim J.Genomics 71:222-234(2001)Research Topic:Cancer

Note: This product is for in vitro research use only