

Human IgG

Catalog No: #IG701



Package Size: #IG701-1 100mg #IG701-2 200mg #IG701-3 1g

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human IgG
Brief Description	Human IgG
Purification	Affinity chromatography or caprylic acid-ammonium sulfate,95%
Concentration	mg/ml
Storage	Store at -20Λ C

Product Description

Source: Normal, healthy, unimmunized rabbit serum

Appearance: Liquid

Filtration: 0.2um

Buffer :10mM Phosphate buffered salinePBSpH7.4

Preservative: None

Stability: Stable for 24 months

Precautions: This product can be kept at 2-8Λ C for maximum 30 days. Avoid repeated freezing and thawing. If multiple usages are needed, it is recommended to be stored separately in smaller volumes.

Note: This product is for in vitro research use only