

FishB ThyrotropinB ReleasingB HormoneB (TRH)B ELISAB Kit



Catalog No: #EK20148

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #EK20148-1 48T #EK20148-2 96T

Description

Product Name	FishB ThyrotropinB ReleasingB HormoneB (TRH)B ELISAB Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Fish

Application Details

Detect Range:0.5 - 100 pg/ml
Sensitivity:0.4pg/ml
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:competitive inhibition

Note: This product is for in vitro research use only