

Mouse Phospho-extracellular signal-regulated kinase (PERK) ELISA Kit



Catalog No: #EK20227

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20227-1 48T #EK20227-2 96T

Description

Product Name	Mouse Phospho-extracellular signal-regulated kinase (PERK) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:31.25-2000pg/ml
Sensitivity:9.1pg/ml
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only