

Rat Asporin (ASPN) ELISA Kit

Catalog No: #EK20251



Package Size: #EK20251-1 48T #EK20251-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Asporin (ASPN) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rt

Application Details

Detect Range:0.156-10ng/ml

Sensitivity:0.06ng/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only