

Mouse Pleckstrin homology-like domain family A member 2 ELISA Kit

Catalog No: #EK20151

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20151-1 48T #EK20151-2 96T

Description

Product Name	Mouse Pleckstrin homology-like domain family A member 2 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Ms

Application Details

Detect Range:31.2-2000pg/ml
Sensitivity:11.2pg/ml
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only