

# Recombinant 2019-nCoV S Protein RBD (K417N, E484K, N501Y, Mammalian, C-His)


Catalog No: #AG0071

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|               |   |
|---------------|---|
| Product Name  | Recombinant 2019-nCoV S Protein RBD (K417N, E484K, N501Y, Mammalian, C-His)   |
| Purification  | >90% as determined by SDS-PAGE  |
| Other Names   | RBD Protein, Spike RBD Protein, South African (ZA) variant (B.1.351), Beta  |
| Calculated MW | 35kDa   |
| Formulation   | Supplied as solution form in PBS, pH7.5 or lyophilized from PBS, pH7.5  |
| Storage       | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt. |

## Images



## Product Description

In general, proteins are shipped out with blue ice unless customers require otherwise.

Note: This product is for in vitro research use only