

Recombinant 2019-nCoV S Protein RBD (K417T, E485K, N501Y, Mammalian, C-His)


Catalog No: #AG0072

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant 2019-nCoV S Protein RBD (K417T, E485K, N501Y, Mammalian, C-His)
Purification	>90% as determined by SDS-PAGE
Other Names	RBD Protein, Spike RBD Protein, Brazil variant (P.1), Gamma
Calculated MW	35kDa
Formulation	Supplied as solution form in PBS, pH7.5 or lyophilized from PBS, pH7.5
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.

Images



Product Description

In general, proteins are shipped out with blue ice unless customers require otherwise.

Note: This product is for in vitro research use only