

Recombinant 2019-nCoV S Protein RBD (K417N,L452R, T478K, Mammalian, C-His)


Catalog No: #AG0074

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|---------------|--|
| Product Name | Recombinant 2019-nCoV S Protein RBD (K417N,L452R, T478K, Mammalian, C-His) |
| Purification | >90% as determined by SDS-PAGE |
| Other Names | RBD Protein, Spike RBD Protein, IN variant (AY.1), Delta Plus |
| Calculated MW | 35kDa |
| Formulation | Supplied as solution form in PBS, pH7.5 or lyophilized from PBS, pH7.5 |
| Storage | Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt. |

Images



Product Description

In general, proteins are shipped out with blue ice unless customers require otherwise.

Note: This product is for in vitro research use only