

# Recombinant 2020-nCoV S Protein RBD (L452Q, F490S, Mammalian, C-His)


Catalog No: #AG0077

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant 2020-nCoV S Protein RBD (L452Q, F490S, Mammalian, C-His)
Brief Description	Recombinant 2020-nCoV S Protein RBD (L452Q, F490S, Mammalian, C-His)
Purification	>90% as determined by SDS-PAGE
Other Names	RBD Protein, Spike RBD Protein, C.37, Lambda
Calculated MW	35kDa
Formulation	Supplied as solution form in PBS, pH7.5 or lyophilized from PBS, pH7.5
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week . Store at -20 to -80 °C for twelve months from the date of receipt.

## Images



## Product Description

shipping: In general, proteins are shipped out with blue ice unless customers require otherwise.

Note: This product is for in vitro research use only