

Purified anti-human CD16

Catalog No: #28016




Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

Product Name	Purified anti-human CD16
Clone No.	4ACB16
Isotype	Mouse IgG1
Applications	FC
Species Reactivity	Human
Concentration	1mg/ml
Formulation	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide.
Storage	Store at 4°C.DO NOT FREEZE.LIGHT SENSITIVE MATERIAL.

Images



Note: This product is for in vitro research use only