

Purified anti-human CD56

Catalog No: #28017


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Purified anti-human CD56
Clone No.	4AMY31
Isotype	Mouse IgG1,k
Applications	FC
Species Reactivity	Human
Formulation	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide.
Storage	Store at 4°C.DO NOT FREEZE.LIGHT SENSITIVE MATERIAL.

Images



Note: This product is for in vitro research use only