

MouseB 8-isoB prostaglandinB F2 α B (8-iso-PGF2 α)B ELISAB Kit



Catalog No: #EK20158

Package Size: #EK20158-1 48T #EK20158-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	MouseB 8-isoB prostaglandinB F2 α B (8-iso-PGF2 α)B ELISAB Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:31.2-2000pg/ml

Sensitivity:12.4pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only