

Recombinant human Methionine synthase reductase

Catalog No: #AP71604



Package Size: #AP71604-1 20ug #AP71604-2 100ug #AP71604-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant human Methionine synthase reductase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-725aaSequence Info:Full Length
Accession No.	Q9UBK8
Uniprot	Q9UBK8
GeneID	4552;
Calculated MW	107.4 kDa
Tag Info	N-terminal GST-tagged
Target Sequence	MGAASVRAGARLVEVALCSFTVTCTCLEVMRRFLLLYATQQGQAKAIAEEICEQAVVHGFSADLHCISESDKYDL KTETAPLVVVVSTTGTGDPPTARKFVKEIQNQTLPVDFFAHLRYGLLGLGDSEYTYFCNGGKIIDKRLQELGA RHFYDTGHADDCVLELVVEPWIAGLWPALRKHFSSRGQEEISGALPVASPASSRTDLVKSELLHIESQVEL LRFDDSGRKDSEVLKQNAVNSNQSNVIEDFESSLTRSVPLSQASLNIPGLPPEYLQVHLQESLQEEESQVS VTSADPVFQVPISKAVQLTTNDAIKTTLLVELDISNTDFSYQPGDAFVSVCPSNSDSEVQSLQLQLEDKREHCV LLKIKADTKKKGATLPQHIPAGCSLQFIFTWCLEIRAIPKKAFLRALVDYTSDSAERKRLQELCSKQGAADYSRF VRDACALLDLLLAFSPCQPPLSLLLEHLPKLQPRPYSCASSLFHPGKLFHFVFNIVEFLSTATTEVLRKGVCTG WLALLVASVLQPNHASHEDSGKALAPKISISPRTTNSFHLPDDPSIIMVGPPTGIAPFIGFLQHREKLQEHP DGNFGAMWLFPGCRHKDRDYLFKELRHLKHGILHLKVSFSRDAPVGEELAPAKYVQDNIQLHGQQVARIL LQENGIHYVCGDAKNMAKDVHDALVQISKEVGVKLEAMKTLATLKEEKRYLQDIWS
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Involved in the reductive regeneration of cob(I)alamin (vitamin B12) cofactor required for the maintenance of methionine synthase in a functional state. Necessary for utilization of methyl groups from the folate cycle, thereby affecting transgenerational epigenetic inheritance. Folate pathway donates methyl groups necessary for cellular methylation and affects different pathways such as DNA methylation, possibly explaining the transgenerational epigenetic inheritance effects.

References

Restricted role for methionine synthase reductase defined by subcellular localization. Froese D.S., Wu X., Zhang J., Dumas R., Schoel W.M., Amrein M., Gravel R.A. Mol. Genet. Metab. 94:68-77(2008) Research Topic: Metabolism

Note: This product is for in vitro research use only