

Recombinant human Hematological and neurological expressed 1 protein



Catalog No: #AP71605

Orders: order@signalwayantibody.com

Package Size: #AP71605-1 20ug #AP71605-2 100ug #AP71605-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant human Hematological and neurological expressed 1 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:2-154aaSequence Info:Full Length
Other Names	Androgen-regulated protein 2
Accession No.	Q9UK76
Uniprot	Q9UK76
GeneID	51155;
Calculated MW	31.9 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	TTTTTFKGVDPNSRNSRVLRLPPGGGSNFSGLGFDEPTEQPVRKNKMASNIFGTPEENQASWAKSAGAKSSG GREDLESSGLQRNSSEASSGDFLDLKGEGDIHENVDTDLPGLGQSEEKPVPAAPVPSVPAPVPSRRNP PGGKSSLVLG
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

Cloning, expression and subcellular localization of HN1 and HN1L genes, as well as characterization of their orthologs, defining an evolutionarily conserved gene family.Zhou G., Wang J., Zhang Y., Zhong C., Ni J., Wang L., Guo J., Zhang K., Yu L., Zhao S.Gene 331:115-123(2004)Research Topic:Epigenetics and Nuclear Signaling

Note: This product is for in vitro research use only